Description

Evetric is a light-curing, radiopaque nano-hybrid composite for restorative treatment. Evetric cures with light in the wavelength range of 400–500 nm.

Composition

Evetric consists of dimethacrylates (19-20 wt.%). The fillers contain barium glass, ytterbium trifluoride, mixed oxide and copolymers (80-81 wt.%). Additives, initiators, stabilizers and pigments are additional contents (< 1 wt.%). The total content of inorganic fillers is 55–57 vol.%. The particle size of inorganic fillers is between 40 nm and 3000 nm.

Indications

- Class I-V restorations
- Restoration of deciduous teeth
- Extended fissure sealing
- Direct veneers
- Splinting of mobile teeth
- Blocking out of undercuts
- Repair of composite and ceramic veneers

Contraindications

The material should not be used

- If a dry working field cannot be established or the stipulated application technique cannot be applied.
- If the patient is known to be allergic to any of the ingredients.

Side effects

In individual cases, components of Evetric may lead to sensitization. Evetric should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp and cover with a suitable cavity liner).

Interactions

- Materials containing eugenol or clove oil may inhibit the polymerization of composite materials. Consequently, the application of such materials together with Evetric must be avoided.
- Discoloration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Application**1. Shade determination**

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

2. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

3. Cavity preparation

Prepare the cavity according to the requirements of the adhesive technique.

Subsequently, remove all residue in the cavity with water spray and dry with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Only cover very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently use a pressure-resistant cement (e.g. a glass ionomer cement, such as Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Matrix / Interdental wedge

Place a matrix and an interdental wedge, if required.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instruction for Use of the product in use. Ivoclar Vivadent recommends using the self-etching adhesive Evetric® Bond. The bond to enamel can be improved by selectively etching the enamel. Unprepared enamel surfaces must be conditioned with phosphoric acid (e.g. Total Etch).

7. Application of Evetric

- Apply Evetric in layers of max. 2mm or 1.5mm (Dentin shade) and contour/adapt the material to the cavity walls using a suitable instrument (e.g. OptraSculpt®).
 - Sufficient exposure to the curing light prevents incomplete polymerization. For the recommendations regarding exposure time and light intensity see table 1.
 - When using a metal matrix, additionally polymerize the composite material from the buccal or lingual/palatal aspect after removing the matrix, if no Bluephase polymerization light is used or the light probe cannot be ideally positioned, e.g. distant to composite or diverging scattering angle.
 - A flowable composite can be used as a thin initial layer. This optional layer has to be cured separately (please refer to the respective Instruction for Use).
- 8. Finishing / Checking the occlusion / Polishing**
- Finish the restoration with suitable finishers or fine diamonds. Subsequently, check the occlusion. Polish with silicone polishers (e.g. OptraPol®).

Additional information

1. If fresh Evetric is applied to polymerized material, it must first be roughened with diamond grinders, wetted with Evetric Bond and subsequently polymerized.
2. Evetric has to be applied at room temperature. Cold material may be difficult to dispense.
3. Do not use oxidizing disinfectants to disinfect syringes.

Warning

- Avoid contact of unpolymerized material with the skin, mucous membrane and eyes.
- Unpolymerized material may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Storage temperature: 2-28 °C
- Close syringes immediately after usage. Exposure to light causes premature polymerization.

Evetric®

EN Instructions for Use

Light-curing resin-based dental restorative material

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Material kompozit restaurues fotopolimerizues

- Expiry date: see note on syringes and packaging.
- Do not use Evetric after the expiry date.

**Keep out of the reach of children.
For use in dentistry only.**

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

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www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai

Trading Co., Ltd., 2/F Building 1, 881 Wuding Road, Jing An District

200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 3399, Fax +57 1 633 1663

www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118, F-74410 Saint-Jorioz, France

Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, D-73479 Ellwangen, Jagst, Germany

Tel. +49 7961 889 0, Fax +49 7961 6326

www.ivoclarvivadent.de

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate

Veera Desai Road, Andheri (West), Mumbai, 400 053, India

Tel. +91 22 2673 0302, Fax +91 22 2673 0301

www.ivoclarvivadent.in

Ivoclar Vivadent Marketing Ltd.

The Icon, Horizon Broadway BSD, Block M5 No. 1

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www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale, PO Box 303011 North Harbour

Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

ul. Jana Pawla II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

www.ivoclarvivadent.pl

Ivoclar Vivadent LLC

Prospekt Andropova 18 korp. 6/, office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204

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www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Compass Building, Feldspar Close, Warrens Business Park

Enderby, Leicester LE19 4SD, United Kingdom

Tel. +44 116 284 7880, Fax +44 116 284 7881

www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

www.ivoclarvivadent.us

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Ivoclar Vivadent AB