Multilink® Hybrid Abutment

EN Instructions for Use

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

Hybrid abutment crown

Description

Multilink® Hybrid Abutment is a self-curing luting composite for the permanent cementation of ceramic structures made of lithium disilicate glass-ceramic (LS₂) or zirconium oxide on titanium/titanium alloy or zirconium oxide bases (e.g. abutment or adhesive basis) in the fabrication of hybrid abutments or hybrid abutment crowns.

Shades

- HO 0 (High Opacity)
- MO 0 (Medium Opacity)

Working time

The working and curing times depend on the ambient temperature. Once Multilink Hybrid Abutment has been dispensed from the automatic syringe, the following times apply:

| | At room temperature 23 °C ± 1 °C (73 °F ± 2 °F) |
|--------------------------------------|--|
| Working time | approx. 2 minutes |
| Curing time (including working time) | approx. 7 minutes |

The working and curing times become shorter at higher temperatures, while they become longer at lower temperatures.

Mixing ratio

Multilink Hybrid Abutment is extruded from the automatic syringe in an optimum 1:1 mixing ratio.

Composition

The monomer matrix is composed of dimethacrylate and HEMA. The inorganic fillers are composed of barium glass, ytterbium trifluoride, spheroidal mixed oxide and titanium oxide. The particle size ranges between 0.15 and 3.0 µm. The average particle size is 0.9 µm. The total volume content of inorganic fillers is approx. 36 %.

Indication

Extraoral, permanent cementation of ceramic structures made of lithium disilicate glass-ceramic (LS₂) or zirconium oxide on titanium/titanium alloy or zirconium oxide bases.

Rx ONLY

Date information prepared:

2013-12-20/Rev. 2
649679/WE3

CE 0123

 **Manufacturer**
Ivoclar Vivadent AG
FL-9494 Schaan/Liechtenstein
www.ivoclarvivadent.com

ivoclar
vivadent®
clinical

Contraindication

The product is contraindicated

- for intraoral luting;
- if the stipulated working technique cannot be applied;
- if the patient is known to be allergic to components of Multilink Hybrid Abutment.

Side effects

Systemic side effects are not known. In rare cases allergic reactions to individual components have been reported.

Interactions

Phenolic substances (e.g. wintergreen oil, eugenol) inhibit polymerization. Therefore, materials containing these components must not be used. Oxidative disinfectants (e.g. hydrogen peroxide) may interact with the initiator system and impair the setting process. As a result, the automix syringe must not be disinfected with oxidative agents. The syringe can be disinfected by wiping it with medical alcohol for example.

Application

The contact surfaces must be meticulously prepared in order to ensure an optimum adhesive bond between the base and the ceramic structure. The procedure is described below. The same procedure is used for the Hybrid Abutment and the Hybrid Abutment Crown.

Hybrid abutments and hybrid abutment crowns made of IPS e.max lithium disilicate glass-ceramics (LS₂)

| | Ceramic structure (LS ₂) | Base |
|-----------------------------------|--|---|
| Blasting | – | According to the instructions of the manufacturer |
| Etching | The bonding area is etched with PS [®] Ceramic Etching Gel for 20 s | – |
| Conditioning | Bonding surface with Monobond Plus for 60 s | |
| Adhesive cementation | Multilink Hybrid Abutment | |
| Covering the cement joint | Glycerine gel, e.g. Liquid Strip | |
| Curing | Approx. 7 min auto-polymerization | |
| Polishing the cement joint | Conventional polishers for ceramic/composite resin | |

Hybrid abutments made of zirconium oxide (ZrO₂)

| | Ceramic structure (ZrO ₂) | Base |
|-----------------------------------|--|---|
| Blasting | According to the instructions of the manufacturer | According to the instructions of the manufacturer |
| Conditioning | Bonding surface with Monobond Plus for 60 s | |
| Adhesive cementation | Multilink Hybrid Abutment | |
| Covering the cement joint | Glycerine gel, e.g. Liquid Strip | |
| Curing | Approx. 7 min auto-polymerization | |
| Polishing the cement joint | Conventional polishers for ceramic/composite resin | |

1. Preparation of the base

The following procedure should be observed in the preparation of the base for the cementation with the ceramic structure:

- The base should be prepared according to the instructions of the manufacturer.
- The base is cleaned in an ultrasonic bath or with a steam cleaner and then dried with blown air.
- The base is screwed to the model analog.
- The ceramic structure is placed on the base and the relative position is marked with a waterproof pen. This makes it easier to attain the correct relative end position when the parts are assembled at a later stage.
- The emergence profile of the base must not be blasted or modified in any way!
- **If the manufacturer recommends that the contact surfaces be sandblasted, the following procedure should be observed:**
 - Hard modelling wax is applied to protect the emergence profile, since this type of material is easy to remove later on.
 - The screw channel is also sealed with wax.
 - The bonding area surface is carefully blasted according to the instructions of the manufacturer.
 - An instrument and a steam jet are used for cleaning. The wax has to be removed meticulously.
- After the bonding surface has been cleaned, it must not be contaminated under any circumstances, as this would impair the bond.
- Monobond Plus is applied on the clean bonding surface and allowed to react for 60 s. After the reaction time, any residue is dispersed with air that is free of water and oil.
- The screw channel is sealed with a foam pellet or wax. The bonding surface must not be contaminated in the process.

2. Preparation of the ceramic structure

2.1. IPS e.max lithium disilicate glass-ceramic (LS₂)

The following procedure must be followed in the preparation of the ceramic structure for the cementation with the base:

- The ceramic structure **must not** be blasted when it is prepared for cementation.

- The ceramic structure is cleaned in an ultrasonic bath or with a steam cleaner and then dried with blown air.
- After the cleaning step, the bonding surface must not be contaminated, since this would negatively influence the bond.
- Wax can be applied to protect the outer surfaces and the glazed areas.
- The bonding surface is etched with 5% hydrofluoric acid (IPS Ceramic Etching Gel) for 20 s.
- Subsequently, the bonding surface is thoroughly rinsed under running water and then dried with oil-free air.
- Monobond Plus is applied to the clean bonding surface and left to react for 60 s. After the reaction time, any residue is distributed with blown air that is free of water and oil.
- Excess in the screw channel is carefully removed with, e.g. a Microbrush or brush using rotary movements.
- **Important: Excess must not be removed before curing has started, i.e. 2–3 minutes after mixing. For the purpose, a suitable dental lab instrument (e.g. Le Cron) is used and the components are held in place with light pressure in the process.**
- Glycerine gel is applied (e.g. Liquid Strip) on the cementation joint to prevent the formation of an inhibition layer. The glycerine gel must be applied cautiously to avoid blending it with or displacing the composite. The gel must be left on the cementation joint until polymerization is complete.

2.2 Preparing ceramic structures made of zirconium oxide (ZrO₂)

The following procedure must be observed in the preparation of the ceramic structure for the cementation with the base:

- In order to facilitate its handling, the ceramic structure is attached to a holder (e.g. Microbrush) before it is blasted.
- The outer surfaces of the ceramic structure are protected with a wax coating.
- In order to better control the blasting procedure, the bonding surface can be coloured with a permanent marker.
- The bonding surface of the ceramic structure is cleaned by means of blasting with Al₂O₃ 100 µm at max. 1 bar and at a distance of max. 1–2 cm from the blasting nozzle (**The instructions of the manufacturer must be observed**).
- The ceramic structure is cleaned with a steam jet or in an ultrasonic bath.
- After cleaning, the bonding surface must not be contaminated under any circumstances, since this may negatively influence the bonding results.
- Monobond Plus is applied to the clean bonding surface and allowed to react for 60 s. After this reaction time, any residue is dried with air that is free of water and oil.

3. Cementation with Multilink Hybrid Abutment

The following instructions must be observed in the cementation procedure:

- The cleaned and conditioned components (ceramic structure, base) are laid out ready for cementation.
- **The subsequent cementation procedure must be carried out quickly and without interruption. The working time of Multilink Hybrid Abutment is approx. 2 min. at 23 °C/73 °F (± 1°C/2 °F).**
- As a general rule, a new mixing tip is attached to the Multilink Hybrid Abutment syringe prior to each use.
- A thin layer of Multilink Hybrid Abutment is directly applied from the mixing syringe to the bonding surface of the base and **the bonding surface of the ceramic structure**.
- The mixing tip is left on the Multilink Hybrid Abutment syringe until the next use. The remaining cement polymerizes in the tip and functions as a seal.
- The ceramic structure is placed on the base in such a way that the position markings are aligned.
- The parts are lightly and evenly pressed together and the correct relative position of the components is checked (transition between the base and ceramic structure).
- Subsequently, the parts are tightly pressed together for 5 s.

- Next, the luting composite is auto-polymerized within 7 min.
- **Important: The parts should not be moved until Multilink Hybrid Abutment has completely cured. They can be held immobile with e.g. diamond-coated tweezers.**
- After the completion of auto-polymerization, the glycerine gel is rinsed off with water.
- **The cementation joint should be cautiously polished with rubber polishers at a low speed (< 5,000 rpm) to avoid overheating.**
- Any cement residue left in the screw channel is removed with suitable rotating instruments.
- The restoration is steam cleaned.

Additional notes

- Multilink Hybrid Abutment should be processed at room temperature. Low(er) temperatures may cause difficulty in dispensing and mixing the material and may prolong the working and curing time.
- Unpolymerized material may stain clothing.

Warning

Contact of uncured Multilink Hybrid Abutment with skin/mucous membrane and eyes must be prevented. In an unpolymerized state, Multilink Hybrid Abutment may cause slight irritation and lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage


- Multilink Hybrid Abutment must not be used after the indicated date of expiration.
- Storage temperature of Multilink Hybrid Abutment: 2–28 °C/36–82 °F
- Shelf life: see note on the syringe
- The mixing tip is left on the syringe after each use. The material left in the tip will cure and the mixing tip will therefore serve as a seal.

Keep out of the reach of children!


For use in dentistry only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.