LEDition®

**Operating Instructions**

Page 2

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Seite 5

Mode d'emploi

Page 9

Istruzioni d'uso

Pagina 11

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Pagina 14

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Página 17

Bruksanvisning

Sida 20

Brugsanvisning

Side 23

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Sivu 26

Bruksanvisning

Side 29

Bedieningshandleiding

Pagina 32

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Page / Seite 50

For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist

Made in Austria

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

ivoclar[®]
vivadent[®]
clinical

Safety

Intended use

The LED polymerization light produces energy-rich blue light. It is intended for the polymerization of light-curing dental materials directly at the dental unit. These Operating Instructions provide recommendations and guidelines for the intended use of the light. Please read them carefully.

Indication

The LED polymerization light is particularly suitable for the polymerization of light-curing dental materials activated in the wavelength range of 430–490 nm.

Contraindication

Materials, the polymerization of which is activated outside the wavelength range of 430–490 nm. If you are not sure about certain products, please ask the manufacturer of the corresponding material.

This unit should not be used near a flammable anaesthetic mixture with air, oxygen or nitrous oxide.

Portable and mobile HF communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.

Patients with heart or brain pacemakers can be treated with LEDition. Investigations have shown that there are no (electromagnetic) incompatibilities between the curing light and pacemakers.

Signs and symbols

Instructions for Use

Contraindication

Symbols on the curing light

Double insulation (Apparatus complies with safety class II)

Protection against electrical shock (BF type apparatus)

Use only in protected areas

Observe Operating Instructions

Do not dispose of the apparatus with normal domestic waste. Please visit your local Ivoclar Vivadent website for detailed information about the disposal of the unit.

NRTL
(Nationally Recognized
Testing Laboratory)

Safety notes

The apparatus has been designed according to IEC60601-1, EN60601-1 and UL60601-1; CSA C22.2 No. 601.1 (2003), complies with the relevant EU regulations, and has been certified by TÜV Product Service as an internationally accredited testing body. The apparatus has been shipped from

the manufacturer in a safe and technically sound condition. In order to maintain it in this condition and to ensure riskfree operation, the notes and guidelines in these Operating Instructions have to be observed. To prevent damage to equipment and protect patients, users and third parties from risks, the following safety instructions have to be observed:

1. Usage and liability

The LED polymerization light must be employed solely for the intended use. Any other uses are contraindicated. Liability cannot be accepted for damage resulting from misuse or failure to observe the Operating Instructions. Furthermore, the user is responsible for testing the LED polymerization light for its use and suitability for the intended purposes. This is particularly important if other equipment is utilized in the immediate vicinity of the apparatus at the same time as the curing light is being used. Use only original components. The manufacturer does not accept any liability for damage resulting from the use of other components.

2. Operating voltage

Before switching on the unit, make sure that

- a) the local power supply complies with the voltage indicated on the rating plate and
- b) the unit has acquired ambient temperature.

3. Indications of compromised safety

If there is any doubt about safe operation being certain, the unit must be disconnected from the power supply and secured against accidental operation. Safe operation may be impaired if, for example, the apparatus is visibly damaged or no longer works correctly. The unit is only completely disconnected from the power supply if the power cord has been pulled out of the plug-in socket.

4. Disinfection and sterilization

Contaminated surfaces of the apparatus, light probes and anti-glare cone must be disinfected or autoclaved before each use.

5. Light output

Before each use, make sure that the light output permits adequate curing. For this purpose, check the light probe for contamination and damage.

6. Heat development

As it is the case with all polymerization lights, the high light output results in a certain development of heat. Prolonged exposure of the pulp and soft tissues may result in damage. Therefore, uninterrupted curing times of more than 40 seconds in the same area, as well as direct contact with the gingiva, oral mucous membrane, or skin have to be prevented. Prevent heat accumulation when polymerizing with a rubber dam in place.

7. Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant for the eyes and may result in injury. It is therefore advisable to use the anti-glare cones provided. Individuals who are generally sensitive to light, who take photosensitizing drugs, or who have undergone surgery and people who work with the apparatus or in its vicinity for long periods of time should not be exposed to the light of this device and wear protective goggles that absorb light below 515 nm.

8. Maintenance and cleaning

Any repair work must only be performed by a certified service center. Liquids or any foreign substances must not enter the handpiece, handpiece holder and particularly the power pack during cleaning (electrical shock hazard).

9. Disposal

Do not dispose of the apparatus with normal domestic waste. Dispose of curing lights according to the corresponding national rules and regulations.

Start-up

Start

The light is activated by pressing the start button and switches off automatically after 2 minutes. If you wish to switch off the light before the set curing time has elapsed, press the start button again. The fan is switched on at the same time as the light is activated. After the curing cycle has ended, the fan continues to run for some more time. The power pack must not be disconnected from the power supply as long as the fan is running.

Curing times

Observe the Instructions for Use of the material applied when selecting the curing time. In general, a maximum curing time of 20 seconds is valid for light-curing adhesives and a maximum of 30 seconds for composites with a layer thickness of up to 2 mm.

Acoustic signals

Start (Stop)	1 beep
After 10 seconds	1 long beep
After 20 and 80 seconds	1 short beep
After 40 and 100 seconds	2 short beeps
After 60 and 120 seconds	3 short beeps

If preferred, the acoustic signals can be turned off. To turn off the acoustic signals, press and hold the start button for approximately 30 seconds while the light is off. To turn the acoustic signals back on, press and hold the start button for another 30 seconds while the light is off.

Product specifications

Cleaning of the housing

Wipe the handpiece and handpiece holder with a standard aldehyde-free disinfecting solution. Do not clean with highly aggressive solutions (e.g. solutions containing orange oil or solutions with an ethanol content of more than 40%), solvents (e.g. acetone) or pointed instruments, which may damage or scratch the plastic. Clean dirty plastic parts with a soap solution.

Cleaning of the light probe

Clean light probes with a disinfectant and a soft cloth. Carefully remove any possible material residue, such as composite material, from the light probe with e.g. fingernails or a plastic spatula. Do not use sharp or pointed instruments, as they may scratch the surface of the light probe and thus reduce light transmission. Light probes and anti-glare cones may be autoclaved (e.g. 134 °C, 3 bar, 5 minutes). Check the light probes for damage. If the light probe is held against a light source and certain segments appear black, then the glass fibres are broken. In this case, the light probe has to be replaced with a new one.

What if ...?

...the light cannot be activated?

a) 5 short beeps

- Handpiece is overheated:
Allow the handpiece to cool down and try again after a while. If the error reappears, please contact your dealer.

b) 8x 5 short beeps

- Power supply defective:
Check your local power supply.
- Electronic defect:
Please contact your dealer.

c) 10x 5 short beeps

- Electronic defect:
Please contact your dealer.

Delivery form

- 1 Handpiece holder
- 1 Handpiece with power cord and power pack
- 1 Light probe 10 mm, black
- 3 Anti-glare cones
- 1 Set of Operating Instructions

Technical data

Operating voltage handpiece	7.5 VDC
Power pack	100–240 VAC / 50–60 Hz / max. 0.4 A Output 7.5 VDC / 1.7 A Manufacturer: Friwo
Operating conditions	
Temperature	+10 °C to +40 °C
Relative humidity	30 % to 75 %
Ambient pressure	700 hPa to 1060 hPa
Dimensions of handpiece holder	L = 200 mm, W = 120 mm, H = 70 mm
Weight of handpiece holder	245 g
Light source	3W LED
Wave length range	430 – 490 nm
Light output	typically 600 mW/cm ² [min. 500 / max. 900]
Operation	6 min. on / 4 min. off (intermittent)
Dimensions of the handpiece	L = 220 mm, W = 35 mm, H = 160 mm
Weight of the handpiece	200 g
Warranty period	1 year
Transportation and storage conditions	
Temperature	–20 °C to +70 °C
Relative humidity	10 % to 75 %
Ambient pressure	500 hPa to 1060 hPa

Store the apparatus in a closed, roofed room. Protect the device from severe jarring.

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Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.