# Telio® CS Inlay

# Telio® CS Onlay

English

**Brugsanvisning**  
**Käyttöohjeet**  
**Bruksanvisning**  
**Productinformatie**  
**Οδηγίες Χρήσεως**  
**Kullanım kılavuzu**  
**Инструкция**  
**Instrukcja stosowania**

**Instructions for Use**  
**Gebruuchsinformation**  
**Mode d'emploi**  
**Istruzioni d'uso**  
**Instrucciones de uso**  
**Instruções de Uso**  
**Bruksanvisning**

- Lyspolymeriserende provisorisk fyllningsmateriale
- Valokovetteinen väliaikaistätymateriaali
- Lysherdende odontologisk fyllningsmateriale
- Lichtuithardend provisorisch vulmateriaal
- Φωτοπολυμεριζόμενο, προσωρινό εμφρακτικό υλικό
- Isikla sertleşen geçici restorasyon materyali
- Светоотверждаемый временный пломбирочный материал
- Światłoutwardzalne materiały do czasowego wypełniania ubytków

- Light-curing temporary restorative material
- Lichthärtendes provisorisches Füllungsmaterial
- Matériau pour obturation provisoire photopolymérisable
- Materiale da restauro foto-indente per otturazioni provisorie
- Material de obturación provisional fotopolimerizable
- Material fotopolimerizável para restaurações provisórias
- Ljushärdande temporärt fyllnadsmaterial

For dental use only.  
Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein  
Ivoclar Vivadent AG,  
FL-9494 Schaan/Liechtenstein

Complies with:  
ISO 4049, EN ISO 4049

**CE 0123**

**ivoclar**  
**vivadent**  
*clinical*

63523801010E3EH

## Instructions for Use

### Description:

Telio® CS Inlay and Telio® CS Onlay are light-curing, single-component materials for temporary restorations, generally not requiring the application of a temporary cement. Both materials are available in a transparent (T) and universal (Uni) shade and are supplied in syringes and Cavifils.

### Composition:

#### Telio CS Inlay

The monomer matrix consists of monofunctional and difunctional methacrylates (43.4 wt%). The fillers are highly dispersed silicon dioxide and copolymers (55.5 wt%). Fluoride (1500 ppm), catalysts, stabilizers and pigments (0.6 wt%) are additional ingredients.

#### Telio CS Onlay

The monomer matrix consists of monofunctional and difunctional methacrylates (36.3 wt%). The fillers are highly dispersed silicon dioxide and copolymers (62.6 wt%). Fluoride (1500 ppm), catalysts, stabilizers and pigments (0.6 wt%) are additional ingredients.

### Indication:

#### Telio CS Inlay

- Temporization of deep inlay preparations with parallel walls
- Sealing of implant screw access holes
- Relining of prefabricated, temporary crowns and bridges made of polycarbonate or methacrylates.

#### Telio CS Onlay

Temporization of large, shallow cavities (onlay preparations). Please note that after curing, Telio CS Onlay is less elastic than Telio CS Inlay, which may complicate removal if undercuts are present.

## Contraindication:

Telio CS Inlay and Telio CS Onlay

- are unsuitable as temporary crown and bridge materials;
- should not remain in the oral cavity for more than 6 weeks;
- should not be used in patients with a known allergy to any of its ingredients.

## Side effects:

After longer periods of service, the material may demonstrate discoloration. However, this does not compromise the performance of the material. In rare cases, Telio CS Inlay and Telio CS Onlay may cause a sensitizing reaction in patients with a hypersensitivity to any of the materials' ingredients. In these cases, the materials should no longer be used.

## Interactions:

Substances containing eugenol/clove oil may inhibit the polymerization of Telio CS Inlay and Telio CS Onlay. Therefore, Telio CS Inlay and Telio CS Onlay should not be applied in combination with zinc oxide-eugenol cements.

## Application:

### Preparation appointment

1. Prepare, place base or liner and take impression as usual.
2. Clean cavity with water spray.
3. Isolate prepared teeth; use a matrix (e.g. OptraMatrix) if necessary.
4. Application of Telio CS Desensitizer (optional):
  - a) The dentin surfaces should be clean and dry.
  - b) Apply Telio CS Desensitizer and brush the material into the dentin for 10 seconds using a suitable instrument (brush, application brush).
  - c) Carefully disperse excess material/dry with blown air.
5. Place a sufficient quantity of Telio CS Inlay or Telio CS Onlay into the cavity using a spatula or another suitable instrument (e.g. OptraSculpt). If you use Cavifils, the material can be applied directly into the cavity by means of an application instrument (e.g. Cavifil® Injector).  
Contour the restoration, leaving minimum excess.

6. Increments of up to 4 mm can be cured in 10 seconds with a standard curing light (at approx. 650 mW/cm<sup>2</sup>, e.g. bluephase using the LOP mode).
7. Remove matrix if required.
8. Check occlusion and adjust if required.

## Insertion appointment

1. Insert a suitable instrument (e.g. probe/scaler) into the temporary restoration. Remove Telio CS Inlay or Telio CS Onlay from the cavity.
2. Subsequently, clean the cavity (e.g. with rubber cups and oil- and fluoride-free cleaning paste).
3. Try in and cement the permanent restoration in place (with e.g. Multilink Automix, Variolink II).

## Notes:

-  If Telio CS Inlay / Telio CS Onlay is applied from the Cavifil directly into the oral cavity of the patient, we recommend using the Cavifil only once due to hygienic reasons (prevention of cross-contamination between patients).
- Do not use oxidizing disinfectants to disinfect syringes and Cavifils.
- Telio CS Inlay / Telio CS Onlay remains well in place in standard inlay preparations. If Telio CS Inlay or Telio CS Onlay is applied to undercuts of the proximal region, retention can be improved. In this case, do not place a matrix, but a wooden interdental wedge before applying Telio CS Inlay / Telio CS Onlay.
- Telio CS Inlay / Telio CS Onlay bonds with light-curing base/lining materials because of the similar chemical compositions. If such materials are utilized, isolate the base with glycerine gel (e.g. Liquid Strip). This prevents the base/lining from being accidentally removed together with Telio CS Inlay / Telio CS Onlay.
- Use silicone finishers (e.g. Astropol F) or tungsten carbide finishers for grinding and excess removal. Excess material may also be removed with a scalpel.
- Telio CS Inlay / Telio CS Onlay can be contoured more easily by wetting the instrument with Telio CS desensitizer or with an unfilled bonding agent (e.g. Heliobond).

### Warning:

Unpolymerized material should not come in contact with skin/mucous membrane or eyes. Unpolymerized Telio CS Inlay / Telio CS Onlay may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

### Storage

- Do not use Telio CS Inlay / Telio CS Onlay after the indicated date of expiration.
- Shelf life: see note on Cavifils, syringes or packages.
- Close syringes / Cavifils immediately after usage.
- Exposure to light causes premature polymerization.
- Storage temperature: 2–28 °C / 36–82 °F.

### Keep the material out of the reach of children.

For use in dentistry only.

### Date information prepared:

07/2010; Rev. 0

### Manufacturer

Ivoclar Vivadent AG,  
FL-9494 Schaan/Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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## Ivoclar Vivadent UK Limited
## Ivoclar Vivadent, Inc.