Käyttöohjeet
Bruksanvisning
Productinformatie
Οδηγίες Χρήσεως
Kullanım Kılavuzu
Инструкция по применению
Instrukcja stosowania

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Istruzioni d'uso
Instrucciones de Uso
Instruções de Uso
Bruksanvisning
Bruksanvisning

- Itsekovetteinen tai valinnaisesti valokoveteinen yhdistelmämuovisementti restauraatioiden kiinnitykseen implantitajaksiin.
- Kiemisk herdede kompositsement met mulighet for lyshærding til sementering av restaureringer på implantatabuments.
- Zelfhårdend bevestigingscomposit met de optie van lichtuitharding, voor het bevestigen van restauraties op implantaat-abuments.
- Αυτοπολυμερίζουσα ρητίνωδης κοκία με επιλογή φωτοπολυμερισμού για τη στερέωση αποκαταστάσεων σε κολοβάματα εμφυτεύμάτων.
- Restorasyonların implant abutmanlarına simanlanmasında kullanılan, kendinden polimerize olan, ışık polimerizasyon opsiyonlu yapıştırma kompoziti.
- Самоутверждающий композит с опцией световой полимеризации для фиксации реставраций на абзменках имплантатов.
- Cement adhezijnyj o rodwójnym mechanizmie polimeryzacji, do osadzania uzupelnieni protetycznych na implantach.

- Self-curing luting composite with light-curing option for the cementation of restorations on implant abutments.
- Selbsthärrendes Befestigungscomposit mit optionaler Lichthärdung zur Befestigung von Restaurationen auf Implantatabuments.
- Colle auto-adhésive et autopolymérisante avec option photopolymérisation, pour assembler des restaurations sur piliers Implantaires
- Cemento composito autoindurente con fotoindurimento opzionale per la cementazione di restaura su abutment implantari.
- Cemento composito autopolimerizable con opción de fotopolimerización, para la cementación de restauraciones sobre pilares de implantes.
- Composito de cimentação autopolimerizável, com opção fotopolimerizável, para a cimentação de restaurações e pilares de implantes.
- Självhärdande kompositsement med ljushärdningsmöjlighet för cementering av restaurationer på implantatdistanser.
- Selvhardende cement med valgfri lyshærding til cementeringen af restaureringer på implantatabuments.

Instructions for Use

Description

Multilink® Implant is a self-curing luting composite with light curing option for the adhesive cementation of indirect restorations (crowns and bridges) made of metal, metal-ceramic and all-ceramic on implant abutments. Monobond Plus is to be used as a coupling agent to achieve a bond to precious and non-precious alloys as well as all-ceramic materials.

Shades

- MO 0 (white opaque)
- MO 1
- transparent

Working time

The working and curing times depend on the ambient temperature. Once Multilink Implant has been dispensed from the automix syringe, the following times apply:

	At room temperature 23 °C ± 1 °C	Intra-orally approx. 37 °C ± 1 °C
Working time	90 ± 15 sec	90 ± 15 sec
Curing time (without working time)	300 ± 30 sec	120 ± 30 sec

Mixing ratio

When Multilink Implant is dispensed from the automix syringe, the cement is automatically mixed in the optimum ratio of 1 : 1.

Composition

The monomer matrix is composed of dimethacrylate and HEMA. The inorganic fillers are barium glass, ytterbium trifluoride, spheroid mixed oxide. The particle size is 0.25–3.0µm. The mean particle size is 0.9µm. The total volume of inorganic fillers is approx. 40%.

Indication

Multilink Implant is used for the permanent cementation of indirect restorations on implant abutments where an optimum bond is required. It is suitable for:

- crowns and bridges made of
 - metal and metal-ceramic
 - all-ceramic
- abutments made of
 - oxide-ceramics (e.g. zirconium oxide)
 - metal (e.g. titanium)

CE 0123

For dental use only.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein
Ivoclar Vivadent AG,
FL-9494 Schaan/Liechtenstein

ivoclar
vivadent
clinical

Contraindication

The use of Multilink Implant is contraindicated

- if a dry working field cannot be established or the stipulated working technique cannot be applied
- if a patient is known to be allergic to any of the ingredients of Multilink Implant
- if restorations are to be cemented on more than three implant abutments. (Do not cement on more than three abutments at one time.)
- If restorations are to be cemented on natural tooth structure

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to single components have been reported.

Interactions

Phenolic substances (e.g. eugenol) inhibit polymerization. Consequently, the application of materials which contain these substances should be avoided. Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn can impair the curing process. Thus, the automix syringe must not be disinfected using oxidative agents. For disinfection, e.g. an alcohol-based disinfectant spray or wipe can be used.

Application

1. Removal of the temporary restoration

Remove the temporary including the temporary abutment if required.

2. Placement and control of the permanent abutment

Follow the manufacturer's instructions.

To improve control when adjusting the restoration and to facilitate the subsequent removal of excess luting cement, a retraction cord can be placed.

3. Try-in of the restoration and establishment of a dry field

Next, the shade, fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic restorations before they are permanently cemented as there is a risk of fracture. If necessary, make corrections with fine diamonds at medium speed and with slight pressure. Polish ground surfaces. Relative isolation of the treatment field – preferably with OptraGate, cotton rolls, saliva ejector or absorbent pads – is essential during adhesive luting procedures.

4. Pre-treatment of the abutment surface

- Seal the abutment screw channel using e.g. a temporary resin restorative material (Telio® CS Inlay)
- Rinse with water spray.
- Air dry the abutments.
- Apply a thin layer of Monobond® Plus using a microbrush and allow it to react for 60 seconds. Subsequently, dry with water/oil-free air.

5. Pre-treatment of the restoration surfaces

In general, proceed according to the instructions of the manufacturer.

Or:

5.1 Preparation of the restoration

5.1.1 Metal or metal-supported restorations

- Sandblast the internal surfaces of the restoration (sandblasting parameters according to the recommendations of the manufacturer of the restorative material) until a uniformly mat surface is obtained.
- Clean the restoration in an ultrasonic unit for approximately 1 minute if required.
- Rinse the restoration thoroughly with water spray and dry with oil-free air.
- **IMPORTANT!** In order to achieve an optimum bond, do not clean the metal surfaces with phosphoric acid.

5.1.2 Zirconium oxide (e.g. IPS e.max® ZirCAD) and aluminium oxide-based restorations

- Sandblast the internal surfaces of the restoration (sandblasting parameters according to the recommendations of the manufacturer of the restorative material).
- Clean the restoration in an ultrasonic unit for approximately 1 minute if required.
- Rinse the restoration thoroughly with water spray and dry with oil-free air.
- **IMPORTANT!** In order to achieve an optimum bond, do not clean the zirconium oxide surfaces with phosphoric acid.

5.1.3 Lithium disilicate glass-ceramic restorations (e.g. IPS e.max Press, IPS e.max CAD)

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 20 seconds or according to the recommendations of the manufacturer of the restorative material.
- Rinse the restoration thoroughly with water spray and dry with oil-free air.

5.1.4 Glass-ceramic restorations (e.g. IPS Empress®)

- Etch with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 60 seconds or according to the recommendations of the manufacturer of the restorative material.
- Rinse the restoration thoroughly with water spray and dry with oil-free air.

5.2 Apply Monobond Plus to the pre-treated surfaces with a brush or microbrush and allow it to react for 60 seconds. Subsequently, disperse the excess of Monobond Plus with a strong stream of air.

6. Application of Multilink Implant

For each application attach a new automix tip to the syringe. Dispense Multilink Implant from the automix syringe and apply the desired quantity directly onto the restoration. As the luting material will cure in the used mixing tip, it may serve as a seal for the contents of the syringe until needed once again.

Note

Multilink Implant should be used immediately after it has been dispensed from the automix syringe. The restoration should be inserted quickly.

7. Placement of the restoration and removal of excess cement

a) solely self-curing

Seat the restoration in place and fix/hold. Remove excess material immediately with a microbrush/brush/foam pellet/dental floss or an implant scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins).

b) self-curing with additional light-curing

(quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridges)

Seat the restoration in place and fix/hold. Excess material is light-cured with the polymerization light (approx. 650 mW/cm², e.g. Bluephase, LOP mode) for 2–4 seconds per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of approx. 0–10 mm. Subsequently, excess material can be easily removed with an implant scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins). Subsequently, light-cure all margins again for 20 seconds (approx. 1200 mW/cm², e.g. Bluephase, HIP mode)

Note

As with all other composites, Multilink Implant is subject to oxygen inhibition. This means that the surface layer (approximately 100 µm) does not polymerize during curing, as it comes in contact with atmospheric oxygen. In order to avoid this, the restoration margins should be coated with glycerine gel/air block (e.g. Liquid Strip) following the removal of excess. Rinse away the gel after complete polymerization.

8. Finishing of the restoration

- Remove the retraction cord if required.
- Check occlusion and function and adjust if required.
- Polish restoration margins with polishers (Astropol®, OpraPol® or OpraFine) or discs.

Special note

Multilink Implant should be at room temperature when used. Extruding and mixing of the material can be difficult immediately after removing it from the refrigerator.

Warning

- Unpolymerized Multilink Implant paste may cause slight irritations. Avoid contact with the skin, mucous membranes or eyes. If Multilink Implant comes into contact with the eyes, immediately rinse with copious amounts of water and seek medical advice. If the material comes into contact with the skin, rinse with copious amounts of water.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Multilink Implant after the indicated expiration date.
- Store Multilink Implant in a cool place (2–8 °C/36–46 °F).
- Store automix syringe with mixing tip attached after use.
- Shelf life: see expiration date

Keep material out of the reach of children!

For use in dentistry only!

Date information prepared: 07/2012; Rev. 5

Manufacturer:

Ivoclar Vivadent AG

FL-9494 Schaan/Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

The Material Safety Data Sheet is available online at

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.