



**Make it easy on yourself!**

# Monobond Plus

**One Primer for  
all restorative materials**

- Glass-ceramics
- Noble-metal alloys
- Zirconium-oxide ceramics
- Aluminium-oxide ceramics
- Base-metal alloys
- Composites (indirect restorations)



**ivoclar**  
**vivadent**  
passion vision innovation

# Monobond Plus

## One Primer

For the first time ever, a primer is available which you can use to condition **all types of restoration surfaces**: Monobond Plus.

The innovative combination of three different functional groups – silane methacrylate, phosphoric acid methacrylate and sulfide methacrylate – enables you to establish a **strong and durable bond** to any dental restorative material.

## One reaction time

Monobond Plus requires the same reaction time on all the various material surfaces. This makes it even easier for you to cement restorations:

- Prepare the restoration surface
- Apply Monobond Plus and **allow it to react for 60 seconds**
- Disperse any excess with air

**Done!**

## Easy storage

No special storage conditions are required, as Monobond Plus can be stored **at room temperature**.

Moreover, Monobond Plus helps you streamline your inventory: You only require **one primer** for all the different restorative materials.



### Delivery form

Refill

1 bottle Monobond Plus

5 g



**Monobond Plus makes  
your life a little easier –**

**AS FAR AS THE CEMENTATION OF  
RESTORATIONS IS CONCERNED!**

Descriptions and data constitute no warranty of attributes.

Printed in Liechtenstein  
© Ivoclar Vivadent AG, Schaan / Liechtenstein  
627562/0309/e/BVD

**Ivoclar Vivadent AG**  
Bendererstr. 2  
FL-9494 Schaan  
Fürstentum Liechtenstein  
Tel. +423 / 235 35 35  
Fax +423 / 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**ivoclar**  
**vivadent**<sup>®</sup>  
passion vision innovation