

Silamat® S6

Form and function –
perfectly mixed



- Aesthetic design
- Homogeneous mixing results
- Universal use

Silamat® S6

Form and function – perfectly mixed

The Silamat S6 is a modern universal mixing device for amalgam, glass ionomer cements and other predosed dental materials in capsules and injection capsules.



Form and function

The aesthetic design is not only attractive but also functional. Thanks to its rounded shape and the removable cover the unit can be cleaned easily.



Homogeneous mixing results

The patented and time-tested Silamat® technology with the figure-eight mixing pattern and above-average rotation speed of 4500 rpm help to achieve reliable and homogeneous mixing results.



Universal use

The Silamat S6 can be used to mix the contents of all popular capsules and injection capsules in various shapes and sizes.



Easy handling

Simply load the capsule from above with one swift hand movement, set the mixing time, start, and you're done!

Silamat S6

Art. No. 602 286

TECHNICAL DATA

Rotations per minute
Operating voltage
Dimensions
Warranty

4500 rpm
100-240 VAC
LxWxH = 240x230x170
3 years

Descriptions and data constitute no warranty of attributes.
Printed in Switzerland
© Ivoclar Vivadent AG, Schaan / Liechtenstein
611508/0208/e/RDV

Manufacturer and distributor
Ivoclar Vivadent AG
Bendererstr. 2
FL-9494 Schaan
Principality of Liechtenstein
Tel +423 / 235 35 35
Fax +423 / 235 33 60
www.ivoclarvivadent.com

**ivoclar
vivadent**
passion vision innovation