

Ivostar[®] | Gnathostar[®]

Simply good denture teeth



Well-proven and economical –
for satisfied patients



ivoclar
vivadent[®]
passion vision innovation

Ivostar® | Gnathostar®

The Ivostar/Gnathostar tooth line combines valuable properties to a quality product that meets high requirements. Each stage of life presents a wealth of opportunities. High-quality dentures will not keep people from seizing them and thus from enjoying life to the fullest.

True-to-nature moulds

Ivoclar and Gnathostar teeth come in a wide range of moulds that reflects the tooth shapes found in nature.

- Wide range of moulds
- True-to-nature tooth dimensions
- Authentic layering

Expressive aesthetics

- Wide range of shades
- Harmonious cervical transitions
- Lifelike texture and incisal effects

Manifold indications

- Partial dentures
- Combination dentures
- Hybrid dentures
- Complete dentures

Biocompatible material

Ivostar/Gnathostar is based on the tried-and-tested PMMA material.

- Compatibility with oral tissues
- Resistance to wear
- Resistance to plaque
- Stability of shade

Shades

- 20 Chromascope
- 16 A–D
- 4 Bleach BL

Harmonious shades

Ivostar/Gnathostar are coordinated with the 20 shades of the Chromascope shade guide, which is an ideal instrument for modern combined denture prosthetics. Additionally, the teeth are also available in 16 A–D and 4 Bleach BL shades. These shades provide for a true-to-nature, amazingly life-like appearance of the dentures.

Individualization

The SR Chromasit® Creative or SR Adoro® stains permit characterization of the teeth.

Stratos 200 articulator

The teeth and their positioning must exactly comply with the individual jaw movements of the patient. With the articulator, this information may be transferred to the dentures.

Biofunctional Prosthetic System (BPS®)

BPS dentures are fabricated of coordinated quality products. The Ivostar/Gnathostar tooth line is an important component of this system. Ask your dealer for detailed information about BPS.

Delivery form

Ivostar

15 maxillary anterior tooth moulds
6 mandibular anterior tooth moulds

Gnathostar

5 maxillary posterior tooth moulds
5 mandibular posterior tooth moulds

Physical values

Flexural strength (DIN 13'907)	120 MPa
Modulus of elasticity (DIN 13'907)	3,000 MPa
Vickers hardness HV 0.5/30 (PV-324 internal standard)	190 MPa
Ball hardness KDH 36.5/30 (DIN 53'456)	170 MPa
Bond to heat-curing polymers (DIN 13'914)	130 MPa



CE 0123

Descriptions and data constitute no warranty of attributes.
Printed in Liechtenstein
© Ivoclar Vivadent AG, Schaan/Liechtenstein
563348/en/2014-04

Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
Liechtenstein
Phone +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

ivoclar
vivadent®
passion vision innovation