

Programat®

A STORY OF SUCCESS

PrograBase 3

The multifunctional software for your Programat® furnace



Control | Manage | Protocol | Download

PrograBase 3

The multifunctional software for your Programat® furnace

The latest generation of Programat computer software now offers you even more possibilities.

Control

Control your Programat ceramic, press or sintering furnace directly from your computer. Relevant information, such as program number, firing status and current firing temperature, are indicated on the screen.

Manage

You can easily manage your programs via your computer. The resulting data can then be saved to the furnace and/or the computer. You can change the firing or sintering programs as you wish, create new programs and copy existing data from one furnace to another.

Protocol

Protocol the firing cycles you have conducted by means of PrograBase 3. The protocolled firing cycles can be supplemented by additional data, such as patient information, order number or material information. You can save the protocols, export them to a PDF file or print them, if you wish.

Download

Our homepage www.ivoclarvivadent.com can be accessed with one click with the PrograBase 3 software. Here you can download the latest equipment software, Operating Instructions and firing tables free of charge.

System requirements

- Operating system: Windows XP® (including Microsoft .NET Framework 2.0*), Windows Vista®, Windows 7 or Windows 8/8.1
- Screen resolution: 1024 x 768
- Colours: min. 256
- RAM: min. 256 MB
- CPU: 1 GHz
- Min. 11 MB free hard disk space
- Internet connection

*) If Microsoft .NET Framework 2.0 is not available on your computer, it will be installed together with the software. Free space required: about 300 MB.

Required accessories

USB download cable (included in the delivery form of the furnace)



PrograBase 3 supports the following furnaces:



Programat P300/G2



Programat P500/G2



Programat P700/G2



Programat EP 3000/G2



Programat EP 5000/G2



Programat CS



Programat S1/1600



This product forms a part of our Fixed Prosthetics category. All the products of this category are optimally coordinated with each other.